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A blue folder or book cover is centered in the frame. The text "Smart Health Simplified" is printed in white, sans-serif font in the middle of the folder. The background is a blurred, light-colored surface.

Smart Health Simplified

Technical Specifications

Product Name: Eye Massager
Model: R-G006
Input: 5V = 1A
Battery Voltage: DC 3.7V
Battery Capacity: 1500mAh

Manufacturer: Shenzhen Ruiyi Business Technology Co., Ltd.
Qianhai Complex A201, Qianwan Road 1, Qianhai Shenzhen-Hong Kong
Cooperation Zone, Shenzhen, 518000 P.R.China

US Importer: JOICOM CORPORATION
14129 The Merge Street, Building 3 Unit A, Eastvale, CA 92880

UK Importer: EEON(UK)LIMITED
Kemp House 160 City Road, London, EC1V 2NX, United Kingdom

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UK REP Richers & Vision Ltd.
Lion Court, Farnham Road, Bordon, England, GU35 0NF

EU REP C&E Connection E-Commerce (DE) GmbH
Zum Linnegraben 20, 65933, Frankfurt am Main, Germany
Email: info@ce-connection.de

The Bluetooth® word mark and logos are registered trademarks owned by
Bluetooth SIG, Inc. and any use of such marks by RENPHO is under license.
Other trademarks and trade names are those of their respective owners.

FCC ID: 2A26P-RG006 IC: 27707-RG006

This device complies with part 15 of the FCC
Rules. Operation is subject to the following two
conditions: (1) This device may not cause harmful
interference, and (2) this device must accept any
interference received, including interference that
may cause undesired operation.

Made in China
NA2

RENPHO EYERIS 2

| Relaxation for Your Eyes with Ultimate Comfort

Soothing eye massager on acupoints

2 massage & heat levels

Built-in Bluetooth® speakers

Eyeris 2 - Eye Massager

Smart Health Simplified

Refresh your eyes after a long day and enjoy a relaxing massage around the eyes and temples. Soothe away headaches, eye strains, dryness and puffiness after a day at work. Designed to suit your personal preference, the RENPHO Eyeris 2 has two massage and heating levels. The sleek portable design allows you to carry the product on-the-go to the office or to travel!

Massage to De-Stress

2 levels of massage intensity and vibration mode

Heating Function

Choose between low or high heat at
104°F-113°F (40°C-45°C)

Portable Sleek Design

Rechargeable and foldable into 180°

Immersive Music Experience

Connect via Bluetooth® to play your favorite
tunes through built-in speakers

Long Battery Life

Works up to 240 minutes on a full charge

Includes: Eyeris 2, Type-C Charging Cable, Storage Bag

RENPHO EYERIS 2