Mebak

The Key to a Better Life

Mebak

Name: Sea&Mew Consulting GmbH
Address: Mittenhuber Straße 4, 92318 Neumarkt
Contact Person: Yuan Li
Tel: +4915224685061
Email: Compliance.EU@outlook.com

Name: Sea&Mew Accounting Ltd
Address: Electric Avenue Vision 25, London, Enfield EN3 7GD
Contact Person: Billy Han
Tel: +447399648608
E-mail: info@seamew.net

Mebak N1

Product Size: 230*160*140mm

Battery: 7.4V/2000mAh

Charger specifications: Type-C 5V= 2A

Net Weight: 1.26Kg

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
Thank you for purchasing Mebak. This User Manual contains important information about safety, use, and warranty. Before using Mebak products, please carefully read all instructions and warning signs in this User Manual and on Mebak, and properly keep the Manual for future reference.

Safety Statement

- Do not expose the shawl to sun or fire. Keep it away from the heat source.
- Do not immerse the shawl in any liquid and keep it ventilated.
- For adult use only. Children are not allowed to use or play with the shawl.
- Do not use the shawl in sensitive parts of the human body such as the Adam's apple, mouth, spine, heart, or genitalia.
- Do not use the shawl on the same part of the body for an extended period of time continuously in case of pain or bruising caused by excessive massage. Stop using the shawl immediately and seek medical attention promptly if there is any pain or discomfort.
- Do not place your fingers or toes into the gaps during use to avoid pinch injuries.
- Do not use if the massage area had surgery before or has discomfort symptoms such as skin damage, redness, inflammation, or congestion.
- Do not operate continuously for more than 30 minutes. Please ensure that the device rests for 5 minutes before re-use.
- Do not disassemble or change the shawl in any way. Improper disassembly or reassembly may lead to the risk of electric shock, explosion, or contact with battery chemicals.
- Never leave the charging or working shawl unattended.
- Turn off the MEBAK N1 shawl for storage or travel. Keep out of reach of children.
- To protect the environment, the shawl should be handled in accordance with local laws and regulations.

Note: Please consult a doctor or contact technical support if you are not sure if you can use the shawl.

Packing List

Technical Parameters

Battery Capacity	7.4V / 2000mAh
Mode	2 modes (constant speed and rhythmic)
Hot compress	Two levels (low-level 40℃=104℉, high-level 50℃=122℉)
Intensity	Two levels (low and high)
Input	5V $\overline{\rule{0.5em}{0.4pt}}$ 2A
Net weight	2.78 lb (1.26KG)
Dimensions (unfolded)	43.3*6.3*5.5in (1,100*160*140mm)

Operating Instructions

1. Before initial use, please use the attached USB-C cable with a 5V/2A output adapter to fully charge the shawl (The red light is pulsating during charging, and turns into a solid white light after the shawl is fully charged).

2. Function operation:

① Power button:

Press and hold the power button for more than 2s to turn on the shawl, and the default working time is 10 minutes before turning off;

Press and hold the power button for more than 2s to turn off the shawl.

② Mode button:

When the shawl is powered on, the default mode is constant speed mode;

Short press the mode button to switch massage modes: constant speed mode - rhythmic mode - cycle in sequence.

③ Hot compress button:

When the shawl is powered on, the default status is low-level hot compress;

Short press the hot compress button to switch hot compress levels: low-level - high-level - off - cycle in sequence.

④ Intensity button:

When the shawl is powered on, the default intensity is low-level;

Short press the intensity button to switch intensity levels: low-level - high-level - cycle in sequence.

3. Place the shawl around the neck area on the back, and fasten the buckle at the back. With both hands holding the strap, then the shawl can be turned on and used.

4. After being scheduled for 10 minutes, the shawl can be turned on again for use if necessary.

5. The average running time of the shawl is 1-2 hours depending on the speed and the compress pressure applied during use. The average charging time is approximately 2.5 hours.

6. Charging while using is not allowed.

7. Fully charge the shawl for not being used for a long time.

Warranty

IF WITHIN 12 MONTHS FROM THE DATE OF PURCHASE, THIS PRODUCT FAILS DUE TO A DEFECT IN MATERIAL OR WORKMANSHIP, MEBAK LIMITED WILL REPAIR OR REPLACE THE PRODUCT, OR NECESSARY COMPONENTS, FREE OF CHARGE.

THE WARRANTY EXCLUDES:

- a. Damage caused by accidents, abuse, improper handling, or transportation;
- b. Damage caused by unauthorized disassembly or repair;
- c. Damage not used in accordance with MEBAK instructions;
- d. Damage exceeding the product cost;
- e. Damage caused by customers' incorrect storage or protective conditions;
- f. Failure to provide a warranty card or proof of purchase date.

The above restrictions may not apply to you for some regions and countries do not allow restrictions of damages.

This warranty card guarantees specific legal rights, and other rights may vary by country and region, and by state. The warranty service is achieved through the operational policies and procedures of Mebak and all Mebak's international distributors. If you are a customer, please contact the corresponding Mebak's business entity, distributor, or retailer for whom you directly purchased the product for warranty and return merchandise authorization procedures. This warranty is initiated and executed by Mebak.

FAQ

You may encounter the following problems. Please contact customer service or dealer for technical support in case the shawl cannot be used after troubleshooting. Do not attempt to repair the shawl yourself.

Questions	Possible Reasons	Troubleshooting
Unable to charge/ turn on	Incompatible charger	Use the standard power cord with the 5V/2A charger to charge the shawl
	Insufficient power	Charge the shawl with a compatible charger.
Weak massage intensity	Low-level operation	Adjust the intensity level to higher.
	Excessive weight on the massage head	Adjust the body's pressure level
Stops working during use.	Insufficient power	Charge the shawl with a compatible charger.
	Scheduled operation completes.	Turn on the device again for further use.
	Overheat protection works.	Turn on the device again after cooling down for 5 minutes.
No heating	The heating button is not pressed.	Press the heating button.
Noise	Sound emits normally during operation.	It can be used normally if the noise is not particularly loud or abnormal.

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Warranty Card

Mebak Warranty Card	
12 months warranty for non-human damage from the date of purchase	
Recipient's Name	
Contact Information	
Recipient's Address	
Return Reason	☐ return ☐ replace ☐ repair
Order No.	
Note: The above information is only for the use of Mebak after-sales service, and Mebak promises to protect all customer privacy information. Email: amazon@mebakshop.com	
EC REP Name: Sea&Mew Consulting GmbH Address: Mittenhuber Straße 4, 92318 Neumarkt Contact Person: Yuan Li Tel: +4915224685061 Email: Compliance.EU@outlook.com	UK REP Name: Sea&Mew Accounting Ltd Address: Electric Avenue Vision 25, London, Enfield EN3 7GD Contact Person: Billy Han Tel: +447399648608 E-mail: info@seamew.net

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